

MRI/ CT/ X-Ray Testing Summary

Oasis




- MR Safe
- CT / X-Ray Radiopaque

Oasis +



- MR Safe
- CT / X-Ray Radiopaque

Azure




- MR Safe
- X-Ray Radiopaque
- CT Radiolucent (thin gel top lining visible)

Oasis Elite



- MR Safe
- CT / X-Ray Radiopaque (PU expand foam)
- CT / X-Ray Radiolucent (foams invisible, thin gel top lining visible)

ARDS/ Elite Foam




- MR Safe
- CT / X-Ray Radiolucent (foam invisible, similar to air)

MRI: all Pressure Care Positioners are MR safe.

X-Ray: all Pressure Care Positioners are qualitatively visible under C-Arm fluoroscopy.

CT: all Pressure Care Positioners are qualitatively visible under CT imaging and do not demonstrate scattering artifacts.



Tel +353 (01) 451 1755
Email: info@trulife.com

www.Trulife.com