

Custom Pediatric and Neonatal Flowmeter's are Available upon Request.

Pediatric and Neonatal Flowmeter Flow Selection Charts

Single Window View LPM Incremental Flow Selection Chart

Flow Range	Range Code	Flow Settings as They Appear through Window										
1/10 thru 1.5	A3	1.5	1	.90	.80	.70	.60	.50	.40	.30	.20	.10
1/10 thru 1.5	A4	1.5	1.0	0.9	0.8	0.7	0.6	0.5	0.4	0.3	0.2	0.1
1/50 thru 3	A5	3.00	1.00	0.75	0.50	0.30	0.20	0.12	0.08	0.05	0.03	0.02
1/40 thru 2	B1	2.0	1.5	0.9	0.6	0.4	0.25	0.15	0.09	0.06	0.04	0.025
1/32 thru 4	D1	4	3	2	1 1/2	1	3/4	1/2	1/4	1/8	1/16	1/32
1/32 thru 4	D3	4	3	2	1.5	1.0	0.75	0.50	0.25	0.12	0.06	0.03
1/2 thru 6	F3	6.0	5.0	4.0	3.5	3.0	2.5	2.0	1.5	1.0	0.75	0.50
1/2 thru 8	H3	8	7	6	5	4	3	2.5	2	1.5	1	.5
1/2 thru 15	M3	15	10	8	6	5	4	3	2	1.5	1	.5
1/32 thru 15	M4	15	6	2	1 1/2	1	3/4	1/2	1/4	1/8	1/16	1/32
1/4 thru 15	M8	15	12	8	6	4	3	2	1 1/2	1	1/2	1/4
1/16 thru 15	M9	15	10	8	5	3	2	1	0.50	0.25	0.12	0.06
1/8 thru 16	N1	16	8	6	2.5	2	1.5	1	3/4	1/2	1/4	1/8

Double Window View ML Incremental Flow Selection Chart

Flow Range	Range Code	Flow Settings as They Appear through Window										
25 ML thru 1 L	A2	1	500	250	200	175	150	125	100	75	50	25
		L	ML	ML	ML	ML	ML	ML	ML	ML	ML	ML
25 CC thru 3 L	A6	3	2	1.5	1	500	250	200	150	100	50	25
		L	L	L	L	CC	CC	CC	CC	CC	CC	CC
30 ML thru 4 L	D7	4	3	2	1.5	1	750	500	250	125	60	30
		L	L	L	L	L	ML	ML	ML	ML	ML	ML

Custom flow rates and flow ranges from 1/50th of a liter through, and including 60 LPM, available upon request.